# **BIP Central Venous Catheter**

## Bactiguard Infection Protection

<b>English</b>	2
<b>Dansk</b>	10
<b>Nederlands</b>	18
<b>Suomi</b>	27
<b>Français</b>	35
<b>Deutsch</b>	44
<b>Ελληνικά</b>	53
<b>Italiano</b>	63
<b>Norsk</b>	72
<b>Português</b>	80
<b>Español</b>	89
<b>Svenska</b>	98
<b>Türkçe</b>	106
<b>العربية</b>	114

# Instructions for Use

English

## Single-lumen and multiple-lumen BIP Central Venous Catheters

Read the instruction manual before use

Do not use if packaging is damaged

Sterile using ethylene oxide

Keep away from sunlight

For single use only

Latex free

Do not re-sterilize

Manufacturer traceability number

Keep dry

Temperature limit (Storage temperature 4 °C – 25 °C)

**READ ALL INSTRUCTIONS, WARNINGS AND PRECAUTIONS CAREFULLY BEFORE USE. FAILURE TO DO SO MAY LEAD TO SEVERE INJURY TO THE PATIENT OR PATIENT DEATH.**

### DESCRIPTION OF THE DEVICE

The BIP Central Venous Catheter (BIP CVC) is made from flexible tubing with a blue flex tip and is available in single- or multi-lumen versions. Each catheter is coated with Bactiguard® coating, a noble metal alloy that is applied to the catheter surface, and is supplied in a sterile package in an insertion kit.

The Catheter Body has markings to identify the insertion depth. Every ten centimetres from the distal tip there is a single, double or triple line indicating the respective distance. The fifteenth and twenty-fifth centimetres are marked with the respective numbers. Starting from sixteen centimetres away from the tip, every centimetre is marked with a dot.

The BIP CVC is 100% latex free.

### INDICATIONS AND INTENDED USE

BIP CVCs are designed for use in critical care patients to monitor central venous pressure, sample venous blood and administer drugs and solutions intravenously. Multiple-lumen catheters provide multiple access channels to the central venous circulation through a single insertion site, allowing several functions to be performed simultaneously.

Bactiguard®-coated CVCs have been shown to reduce catheter-related infections in high-risk patients, when compared to uncoated CVCs.<sup>1)</sup>

The Bactiguard® coating on the BIP CVC has been shown to reduce bacterial adhesion to the catheter surface in vitro.

1) Goldsmith et al: Prevention of catheter-related infections by silver coated central venous catheters in oncological patients. *Zbl. Bact.* 283, 215–223 (1995). Bactiguard (under the previous company name Metacot) developed a Bactiguard-coated polyurethane CVC together with Fresenius; Bacti-Guard Safe Seldinger, used in the cited clinical study.

After 30 days exposure to artificial blood the adhesion of *Staphylococcus aureus* to the catheter surface was reduced with over 50% as compared to an uncoated catheter.<sup>2)</sup>

The BIP CVC is for short-term use only.

The BIP CVC is not intended to be used to treat existing infections or as a substitute for a tunnelled catheter in patients requiring long-term therapy.

## PRECAUTIONS

- Rx only.
- Do not use the catheter or any other kit component after the expiry date indicated on the packaging.
- This product is designed for single use only. Do not reuse or re-sterilize the catheter or kit components as this may lead to severe patient injury, including death, due to the risk of undesired changes in the physical properties of the devices and possible cross-contamination.
- Do not use the catheter or components if the packaging is opened or damaged, as the contents may lose sterility.
- Do not alter the catheter, guide-wire or any other kit components during insertion, use or removal.
- To minimise the risk of catheter rupture, the infusion pressure should not exceed 40 PSI (2,76 bar).
- To minimise the pressure generated during flushing procedures, use a 10-ml syringe or larger.
- For short-term use of up to 30 days.
- Central venous catheterisation must be performed by a trained physician who is familiar with anatomical

landmarks, the safe technique and potential complications.

- Used catheters must be properly disposed of as biohazard material and processed according to the facility's protocol.
- Be aware of the risks associated with general and local anaesthesia, surgery and post-op recovery.

## WARNINGS

- Do not place the catheter into, or allow it to remain in, the right atrium or right ventricle. Failure to heed this warning may result in severe injury to the patient or the patient's death.
- Due to the risk of exposure to HIV (Human Immunodeficiency Virus) or other blood-borne pathogens, healthcare workers should routinely take universal precautions relating to blood and bodily fluid precautions in the care of all patients. Aseptic technique must be strictly adhered to during any handling of the device.

## SITUATIONS OF GREATER RISK WHERE EXTRA PRECAUTION IS NEEDED AND THE DEVICE MAY BE CONTRAINDICATED

The situation shall be assessed by a trained physician only. The following situations involve greater risk:

- Where the presence of another device-related infection, bacteraemia or septicaemia is known or suspected.
- Where severe chronic obstructive pulmonary disease exists.
- Where previous irradiation of the prospective insertion site has occurred.

*2) Bacterial adhesion to BIP CVC was evaluated using a modified version of the in vitro test for primary adhesion described by Gabriel et al: Effects of silver on adherence of bacteria to urinary catheters: In vitro studies, Current Microbiology, 30, 17-22, (1995), internal report.*

- Where previous episodes of venous thrombosis or vascular surgical procedures have occurred at the prospective placement site.

**If adverse reactions occur after catheter placement, remove the catheter immediately.**

There are no known contraindications due to the Bactiguard® coating.

## SPECIAL PATIENT POPULATIONS

No studies on Bactiguard®-coated central venous catheters have been conducted in pregnant women, paediatric or neonatal patients. The benefits of using the catheter should be weighed against any possible risk.

## COMPLICATIONS

Potential complications associated with use of central venous catheters are listed below;

- Allergic reaction
- Anaphylactic shock
- Bleeding
- Brachial plexus injury
- Cardiac arrhythmia
- Cardiac tamponade
- Catheter damage
- Catheter embolism
- Catheter occlusion
- Catheter erosion through the skin
- Catheter-related bloodstream infection
- Catheter sepsis
- Catheter tip migration
- Death
- Endocarditis
- Extravasation
- Fibrin formation
- Haematoma
- Haemothorax
- Hydrothorax
- Insertion site erythema

- Myocardial damage
- Nerve damage
- Perforation or laceration of vessels or viscus
- Pneumothorax
- Thoracic duct injury
- Thromboembolism
- Tissue necrosis
- Deep vein thrombosis
- Vessel erosion
- Phlebitis

## INSTRUCTIONS FOR USE

### RECOMMENDED INSERTION PROCEDURE:

**CAUTION:** Adhere to aseptic technique to maintain sterility. Follow basic concepts of infection control, take universal precautions and follow the procedures established in institutional protocols.

1. Check the packaging. Do not use the contents if the package is damaged.
2. Prepare and drape the puncture site as required.

**CAUTION:** To reduce the risk of air embolism during catheter insertion, the patient should be positioned in a slight Trendelenburg position insofar as this is tolerated.

3. Prepare the catheter for insertion by flushing the catheter lumen with sterile solution to prime it and ensure it is not blocked. Clamp the extension lines or attach injection ports to the appropriate lumen extensions. Leave the distal lumen extension uncapped for guide-wire passage.
4. Infiltrate the insertion site with local anaesthetic in accordance with the institutional protocol.

### INSTRUCTIONS FOR INTRODUCER NEEDLE

(instructions for the Y-valve needle are given after the storage and shelf life information)

5. Place the introducer needle on the syringe.
6. Locate the vein, insert the needle, and aspirate. Ensure that a good flow of venous blood is established.
7. Remove the syringe. Note: Pulsatile flow exiting the introducer needle usually indicates an unintentional arterial puncture (remove the syringe and apply pressure to the insertion site in order to prevent the risk of a large haematoma).
8. Retract the guide-wire into the guide-wire dispenser to straighten its J-tip. Gently insert the guide-wire tip through the needle into the vessel. Advance the guide-wire to the required depth.

## PRECAUTIONS:

- Exercise care when inserting or withdrawing the guide-wire.
- Review the patient's history before the catheterisation procedure to check for possible implants that may disturb the process.
- Do not withdraw the guide-wire against the needle bevel to avoid potentially severely damaging, or even cutting off, the guide-wire.
- If resistance to the advancement of the guide-wire is encountered, withdraw the guide-wire and gently try to reintroduce it.
- Advancement of the guide-wire into the right heart chamber can cause arrhythmias, right bundle branch block, and vessel wall- atrial- or ventricular perforation.
- While holding the guide-wire in place, remove the introducer needle. Use the centimetre marks on the guide-wire as a reference to adjust the indwelling

length to achieve the desired depth for the indwelling catheter.

**CAUTION:** Maintain a firm grip on the guide-wire at all times to ensure it does not get drawn into the patient's body.

- If needed, insert the vessel dilator over the guide-wire into the blood vessel to enlarge the puncture site. If further enlargement of the puncture site is needed, use the scalpel. Remove the vessel dilator.

**Warning:** Do not leave the vessel dilator in place as an indwelling catheter, since this may cause vessel wall perforation.

**CAUTION:** The physician must be aware of potential air embolism complications associated with leaving open needles or catheters in central venous puncture sites, or as a consequence of inadvertent disconnection.

**Warning:** Do not cut the catheter to alter length.

9. Thread the tip of the catheter over the guide-wire. Holding the catheter near the skin, advance the catheter into the vein with slight twisting motion. A sufficient length of the guide-wire must remain exposed at the hub end of the catheter to allow a firm grip on the guide-wire.
10. Using the centimetre marks on the catheter as positioning reference points, advance the catheter to the final indwelling position.
11. Hold the catheter at the desired depth and remove the guide-wire.


**CAUTION:** If resistance is encountered when attempting to remove the guide-wire after catheter placement, withdraw the catheter 2 or 3 cm and attempt to remove the guide-wire. Applying undue force during guide-wire withdrawal increases the potential for guide-wire or catheter breakage. If resistance is encountered

again, remove the guide-wire and catheter simultaneously.

12. Verify that the entire guide-wire is intact upon removal.
13. Check the lumen placement by attaching a syringe to each lumen extension and aspirate until a free flow of blood is observed. Connect all lumen extensions to appropriate Luer-lock line(s) as required. Unused ports may be "locked" through the injection cap(s) following standard institutional protocols. Slide clamps are provided on lumen extensions to occlude flow through each lumen during line and injection port changes.

### PRECAUTIONS:

- Use only Luer-lock connectors on catheters and accessories to avoid disconnection problems.
- To avoid damage to lumen extensions due to excessive pressure, each clamp must be opened prior to infusion.
- For total occlusion of flow through each lumen, please ensure that the extension line clamps are pressed completely to the opposite end (As shown in picture below).

14. Secure the catheter position by suturing the integral suture hub to the patient skin. Clamp a suture wing over the catheter to ensure another two suture points. Snap a clamp fastener onto the suture wing to secure the catheter before suturing the suture wing and its cover together on the patient's skin.

15. Verify the catheter tip position using a chest x-ray immediately after placement and before the catheter is used for its intended purpose.

**CAUTION:** The catheter should be positioned so that its distal tip is advanced as far as possible into the superior vena cava and as close as possible to the right atrium without touching it. Verify that the catheter lies parallel to the vessel wall. Physicians must be aware of the possible complications caused by a right atrium perforation if the catheter is advanced too deeply.

**CAUTION:** Do not suture directly on the outside diameter of the catheter, to avoid cutting or damaging the catheter or impeding catheter flow.

**Warning:** When inserting the catheter always leave part of it outside of the patient in order to make space for the suture wing. Without the suture wing the catheter cannot be properly secured.

16. Dress the puncture site in accordance with hospital protocol.

17. It is recommended to note the catheter depth in the patient's records.

**CAUTION:** Complications associated with central venous catheters include air embolism, catheter embolism, cardiac tamponade secondary to vessel wall-atrial- or ventricular perforation, septicaemia, and thrombosis. See Section on COMPLICATIONS.

### RECOMMENDED MAINTENANCE:

**CAUTION:** Due to the risks of exposure to blood-borne pathogens or HIV (Human Immunodeficiency Virus), health care workers should routinely take universal precautions relating to blood and bodily fluids in the care of all patients.

1. Indwelling catheters should be routinely inspected for the desired flow rate, security of dressing, correct catheter position, secure Luer-lock

connections and signs of infection at puncture site. Use the centimetre marks to determine whether the catheter position has changed.

**CAUTION:** A catheter that has slightly changed position and slipped out should never be inserted back into its original position, as this might lead to infection. Fixate at new position or exchange catheter.

**CAUTION:** Only an x-ray examination of catheter placement can ensure that the catheter tip has not entered the heart or no longer lies parallel to the vessel wall. If a change in catheter position is suspected, immediately perform a chest x-ray examination to confirm the catheter tip position and relocate as necessary.

2. Maintain the insertion site with regular meticulous dressing changes using aseptic technique in accordance with institutional protocols.
3. To minimise the pressure generated during flushing procedures, use a 10-ml syringe or larger.

**CAUTION:** Do not let solutions containing acetone or alcohol come into contact with the indwelling catheter since these may weaken the catheter and can cause it to leak or break.

**CAUTION:** The use of a syringe of less than 5 ml could cause the catheter to burst and/or the tip to detach.

4. The physician should evaluate the length of time the catheter is to be left in place to reduce the risk of contamination at femoral access sites.
5. To take blood samples, temporarily shut off port(s) through which solutions are being infused.

**CAUTION:** Compression or kinking of the catheter or extension tubing may cause sudden increases, decreases, or small boluses of vasoactive drugs to be injected

into the patient, which could result in sudden changes in cardiac output.

### **CATHETER REMOVAL:**

1. Remove the dressing.

**CAUTION:** Avoid removing the dressing with scissors or other sharp instruments that may cut the catheter.

2. Remove all sutures that hold the catheter in place. Be careful not to cut the catheter.
3. Slowly remove the catheter by pulling it parallel to the skin. As the catheter exits the site apply pressure with a dressing.

**CAUTION:** To prevent an air embolism after removal of the central venous catheter, cover the wound with a dressing impermeable to air.

4. Inspect the catheter upon removal to make sure that the entire catheter length has been removed, including the tip.

**Warning:** Do not use excessive force when removing the catheter. If withdrawal cannot be easily completed, a chest x-ray should be performed and further consultation requested as necessary.

**CAUTION:** Do not reuse or re-sterilize the catheter.

### **CATHETER EXCHANGE PROCEDURE:**

1. Use the sterile technique.
2. Remove the catheter clamp and sutures.
3. Proceed in accordance with hospital protocol.

**CAUTION:** Do not cut the catheter since this may result in catheter embolism.

### **PRODUCT DISPOSAL:**

- After use  
Used catheters must be properly disposed of as biohazard material and processed according to the facility's protocol.

- Expired, unused product  
Follow the facility's protocol.

## PACKAGING AND STERILITY

The product is supplied sterile if the package is undamaged and unopened. Do not use if the package is opened or damaged. The catheters are for single use only. Do not clean, reuse or re-sterilize a used catheter.

## STORAGE

The BIP CVC kits should be stored, unopened, in their original packaging in a cool dark place. Store at +4 °C – +25 °C.

## SHELF-LIFE

The recommended shelf-life is indicated on each package.


## Y-VALVE NEEDLE

### Instructions for Use:

The Y-Valve Needle is designed to prevent back flow of blood during guide-wire insertion when using the Seldinger technique.

### Guide-wire Insertion Procedure:

1. Attach the syringe to the straight-through connector on the Y-Valve Needle. Locate the vein, insert the needle, and aspirate. Ensure that a good flow of venous blood is established. Keep the syringe attached to the needle.
2. Stretch the flexible J-tip of the guide-wire by pulling it into the guide-wire dispenser. Attach the dispenser tip to the side connector of the Y-Valve Needle, then advance the guide-wire into the vein.

3. After the guide-wire has reached the desired depth, remove the guide-wire dispenser. Then remove the needle while keeping the guide-wire in place.
4. Insert the dilator and catheter as recommended in the instructions for using the introducer needle.

## RAULERSON SYRINGE

### Instructions for Use:

1. After administering a local anesthetic, locate a vein using the thin-wall introducer needle attached to the Raulerson syringe. The vessel may be prelocated with a smaller needle.
2. Using the Guidewire Advancer, straighten the J tip of the guidewire, if used, and advance through the rear end of the syringe plunger.

**WARNING:** To avoid severing or damaging the wire, do not cut the guidewire to alter its length, nor withdraw the guidewire against the needle bevel.

3. Hold the guidewire in place and remove the Raulerson syringe.

**PRECAUTION:** Maintain a firm grip on the guidewire at all times.

4. If making a skin-nick, enlarge the cutaneous puncture site with the cutting edge of the scalpel positioned away from the guidewire.

5. If using the dilator, pass it over the guidewire to enlarge the site as needed.


**WARNING:** To avoid a possible vessel wall perforation, do not leave the dilator in place as an indwelling catheter.

6. Thread the tip of the catheter into a vessel using the guidewire. Holding the catheter near the skin, advance into the vein with a slight twisting motion.

7. Advance the catheter into the final indwelling position. Hold the catheter and remove the guidewire. Check the lumen placement by aspirating through the lumen extensions. Apply dressing per hospital protocol. Verify the catheter tip position by X-ray (or other method in compliance with hospital protocol).

## NEEDLE STOPPER

A needle stopper is included to temporarily shield the sharpened tips of bevelled needles before they are disposed of in accordance with local, state and federal regulations.

**CAUTION:** US Federal Law restricts this device to sale by or on the order of a physician or licensed practitioner.

## PRODUCT INFORMATION

For further information or assistance relating to Bactiguard products, please contact:

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#### **Bactiguard AB**