BIP Endotracheal Tube Evac

Bactiguard Infection Protection

| | |
|-------------------|-----------|
| English | 2 |
| Dansk | 4 |
| Nederlands | 6 |
| Suomi | 8 |
| Français | 10 |
| Deutsch | 12 |
| Ελληνικά | 14 |
| Italiano | 16 |
| Norsk | 18 |
| Português | 20 |
| Español | 22 |
| Svenska | 24 |
| Türkçe | 26 |
| Polski | 28 |
| العربية | 31 |
| References | 33 |
| Symbols | 34 |

Instructions for Use

BIP Endotracheal Tube Evac

Product Description

The BIP Endotracheal Tube Evac (BIP ETT Evac) is a tube designed for insertion through the mouth into the trachea for airway management in medical settings. The BIP Endotracheal Tube Evac has dedicated lumen for subglottic secretion drainage (SSD) and is surface-treated with Bactiguard® coating. The Bactiguard® coating consists of a noble metal alloy which has been shown to significantly reduce bacterial adhesion.

Bacterial adhesion is a contributing factor to bacterial colonization, which in turn increases the risk for Ventilator-Associated Pneumonia (VAP).

The BIP Endotracheal Tube Evac has the following characteristics:

- Three lumen tube made from medical grade polyvinyl chloride (PVC).
- Magill curved with rounded tip and Murphy eye, for oral intubation.
- High volume low pressure cuff and standard connector.
- Radiopaque line and depth marks.
- Separate evacuation lumen which has a dorsal opening above the cuff in subglottic area. Access to the lumen is accomplished via a clear connecting tube with a capped Luer connector.
- Sterile, single use, 100% latex free. Free from phthalates banned by European Union in the Directive 2005/84/EC (DEHP, DBP, BBP, DINP, DIDP, DNOP).

Intended Use

The BIP Endotracheal Tube Evac is indicated for use in airway management by oral intubation of trachea. The evacuation lumen allows for drainage by continuous or intermittent suctioning of contaminated mucous and subglottic secretions that accumulate above the cuff. BIP Endotracheal Tube Evac is intended for patients anticipated to need prolonged mechanical ventilation.

The Bactiguard® coating on the BIP Endotracheal Tube Evac is significantly reducing microbial adhesion. Subglottic secretion drainage (SSD) has been shown effective for prevention of ventilator-associated pneumonia (VAP).

Instructions for Use

The BIP Endotracheal Tube Evac can be used in adults and children above certain size as judged by the physician. Aseptic procedure must be used at all times. For the purpose trained medical staff, following currently accepted medical techniques, should perform intubation and extubation. For the purpose trained medical staff shall determine suitable tube size for each individual patient.

The BIP Endotracheal Tube Evac should be removed/ exchanged after a suitable interval determined individually for each patient, based on clinical indication and according to current guidelines for prevention of infection. If necessary, the BIP Endotracheal Tube Evac can be used for a maximum of 30 days in single or total accumulated use.

Contraindications

The use of a laser beam or electrosurgical active electrode in the immediate area of the device is contraindicated.

Precautions

- Do not use if the packaging is damaged, opened or if the expiration date has passed.
- The BIP Endotracheal Tube Evac is intended for single use only and must not be reused or resterilized. Failure to comply will reduce the product's reliability and functionality and expose patients to risks such as cross infection or injury.
- Non-standard dimensioning of some connectors on ventilators or anaesthesia equipment may make secure mating with the endotracheal tube's 15mm connector difficult.
- Intubation and extubation should be performed using currently accepted medical techniques.
- If stylet is used, pre-test to make sure the stylet can be easily inserted and removed from the endotracheal tube prior to initiation of intubation. The stylet tip may not extend beyond the distal end of the endotracheal tube. If the stylet sheath is torn, cut or lacerated while reshaping or inserting the stylet, do not use for intubation as a damaged sheath can leave residues inside the tube or airways.
- Each tube's cuff, pilot balloon, and valve should be tested by inflation before use. Do not inflate a volume that will deform the cylindrical shape or put excessive pressure on the cuff.
- Avoid damaging the thin-walled cuff during intubation. If the cuff is damaged, the tube must not be used.
- Should extreme chin-to-chest flexing of the head or movement of the patient (e.g. to a lateral or a prone position) be anticipated after intubation, use of a reinforced endotracheal tube should be considered.
- Deflate the cuff prior to reposition of the tube. Movement of the tube with cuff inflated can result in tracheal damage.
- Tubes should be securely anchored to avoid unnecessary tube movement. Perform according to currently accepted medical techniques.
- Seat the connector firmly in both the endotracheal tube and the adapter on the ventilation equipment to prevent disconnection during use.
- A bite block should be used in cases where the patient may bite down and flatten the endotracheal tube.
- Following intubation, inflate the cuff only enough to provide an effective seal at the desired lung pressure. Inflation of cuff by "feel" or by using a measured amount of air is not recommended since resistance is an unreliable guide during inflation. In selecting the sealing pressure, an intracuff pressure measuring device should be used in conjunction with Minimal Occluding Volume or Minimum Leak techniques. Cuff pressure should be monitored. Any deviation from the selected seal pressure should be investigated and corrected immediately.
- Remove syringe from valve housing after cuff inflation. Leaving the syringe will keep the valve open, permitting the cuff to deflate.

- Three-way stopcocks or other devices should not be connected to the inflation valve for extended periods of time. The resulting stress could crack the valve housing and allow the cuff to deflate.
- Expert clinical judgement should be used in preventing, detecting, and correcting possible occlusion of the suction port caused by entrapment of tracheal mucosa. Possible steps to prevent or correct occlusion of the suction port may include monitoring and adjusting cuff pressure according to currently accepted guidelines, periodically releasing suction or switching to intermittent suction if applicable, and using a syringe to administer a bolus of 3–5 cc of air into the suction lumen.
- After moving the patient following intubation, it is essential to verify the correct placement of the endotracheal tube.
- Diffusion of nitrous oxide mixture, oxygen, or air may either increase or decrease cuff volume and pressure. To decrease such diffusion, inflating the cuff with the same gas mixture that will contact the cuff's external surface is recommended.
- The use of Lidocaine Topical Aerosol has been associated with the formation of pinholes in PVC cuffs (Jayasuriya, K.D., and Watson, W.F.: P.V.C. Cuffs and Lignocaine-based Aerosol. *Brit. J. Ann.* 53:1368, 1981). The same authors report that lidocaine hydrochloride solution does not have this effect.
- For lubrication purpose, use only sterile water or water-soluble lubricants. If excessive amounts of lubricants are used, this may partially or totally block the inner lumen of the tube.
- Before extubation, completely deflate the cuff using a Luer-tip syringe until a definitive vacuum is noted.
- After use or in case of expired product, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

ADVERSE REACTIONS

Reported adverse reactions associated with the use of tracheal tubes are many and diverse. Standard textbooks and the scientific literatures should be consulted for specific adverse reaction information.

More detailed discussions of adverse reactions to tracheal tubes are contained in the publications available under Chapter: *References*.

Intubation

1. Remove the BIP Endotracheal Tube Evac from the inner package.
2. Test cuff integrity before intubation. Inflate with a Luer-tip syringe, and then completely evacuate air from the cuff after testing.
3. Ensure that the standard connector is firmly attached to the tube.
4. Intubate the patient using currently accepted medical techniques.
5. Inflate the cuff with sufficient gas mixture to provide an effective seal at the desired lung inflation pressure. Use of Minimal Occluding Volume Technique can reduce occurrence of many of the adverse reactions associated with the use of cuffed tracheal tubes.
6. Remove the Luer-tip syringe from valve housing after cuff inflation. Leaving the syringe will keep the valve open, permitting the cuff to deflate.

Suction

1. Suction of secretion above the cuff using the minimum suction pressure required to effectively remove the secretions. Low continuous suction or moderate intermittent suction may be used. It is suggested that suction levels do not exceed 20 mmHg for continuous suction or 150 mmHg for intermittent suction. The suction lumen should be checked periodically for patency.
2. Follow hospital procedures regarding the suction of the primary airway lumen to remove secretions.

Extubation

1. Before extubation, completely deflate the cuff using a Luer-tip syringe.
2. Extubate using currently accepted medical techniques.

Storage Conditions

- The BIP Endotracheal Tube Evac must be stored in the original packaging at all times.
- Avoid exposure to direct sunlight and ultraviolet light.
- Avoid storage that may deform the device or damage the packaging.
- Store at +4 °C - +25 °C and keep dry.

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References

1. Martin LD, Mhyre JM, Shanks AM, Tremper KK, Kheterpal S. 3,423 Emergency tracheal intubations at a university hospital: airway outcomes and complications. *Anesthesiology* 2011;114(1):42–48.
2. Jaber S, Amraoui J, Lefrant JY, et al. Clinical practice and risk factors for immediate complications of endotracheal intubation in the intensive care unit: a prospective, multiple-center study. *Crit Care Med* 2006;34(9):2355–2361. PMID: 16850003
3. Benjamin B, Holinger LD. Laryngeal complications of endotracheal intubation. *Ann Otol Rhinol Laryngol* 2008;117:1–20.
4. Mort TC. Complications of emergency tracheal intubation: immediate airway-related consequences: part II. *J Intensive Care Med* 2007;22(4):208–215. PMID: 17712056
5. Divatia J, Bhowmick K. Complications of endotracheal intubation and other airway management procedures. *Indian J Anaesth* 2005;49(4):308–318.
6. Haas CF, Eakin RM, Konkle MA, Blank R. Endotracheal tubes: old and new. *Respir Care*. 2014 Jun;59(6):933-52; discussion 952-5. doi: 10.4187/respcare.02868. Review.
7. Efrati S, Deutsch I, Antonelli M, Hockey PM, Rozenblum R, Gurman GM. Ventilator-associated pneumonia: current status and future recommendations. *J Clin Monit Comput*. 2010 Apr;24(2):161-8. doi: 10.1007/s10877-010-9228-2. Epub 2010 Mar 17. Review. Erratum in: *J Clin Monit Comput*. 2012 Dec;26(6):493.
8. Gentile MA, Siobal MS. Are specialized endotracheal tubes and heat-and-moisture exchangers cost-effective in preventing ventilator associated pneumonia? *Respir Care*. 2010 Feb;55(2):184-96; discussion 196-7.

Manufacturer, name and address stated. Anført producent, navn og adresse. Vermelding van fabrikant, naam en adres. Valmistaja, nimi ja osoite ilmoitettu. Fabricant (nom et adresse). Hersteller mit Name und Adresse. Κατασκευαστής, αναφέρεται η επωνυμία και η διεύθυνση. Produttore, con indicazione di nome e indirizzo. Produsent, produsentens navn og adresse. Fabricante, indicação do nome e da morada. Fabricante, denominación y dirección indicadas. Tillverkare, namn och adress. Urećiti, beyan edilen ad ve adres. Producent, nazwa i adres.

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Bactiguard AB